

# Gravity Chair Air-Breathe

model : GC2-001AB

USER'S MANUAL

## USER'S MANUAL



### —Headrest—

An ergonomically wave-shaped foam keeps the head position stable and prevents the head from falling sideways.

### —Seat Liner—

Air-mesh fabric with quilted cushion inside keeps the back dry and comfortable.

### —Wedge Mat—

A Wedge Mat makes sitting angle adjustable to the most comfortable position.

### —PU Wipe-Down Cover—

Fire-retardant and water-repelling treatment makes cover maintenance easy.

### —Seat Depth Spacer (not shown in the photo)—

A foam spacer beneath the cover makes sitting height position adjustable.

### —Side Pocket—

Keeps this User's manual and your small belongings inside.

## ■ GRAVITY SUPPORTS YOU

When you are sitting in the Gravity Chair, four faces, Seat, Back and both Sides of the chair automatically wrap and keep your body posture very stable and comfortable.

By reclining the chair using the Triangle Mat, you will be supported more firmly.

Elbows are pushed a little forward in the chair that automatically makes your shoulders position stay at more natural and relaxed.

## ■ AS YOUR SECOND CHAIR

The Gravity Chair may be your Second Chair for relaxing purpose while the chair you use with upright sitting position for daily dining and studying is your first chair. However, for the people who have difficulty to use ordinary chairs, the Gravity Chair can be their First Chair as it has outstanding capability of holding human body.

## ■ RECLINING

By using Wedge Mat, you can adjust the chair angle to the most comfortable position for each purpose, watching TV set and / or taking a rest. Even with reclining the chair, the angle between seat and back remains right-angled so that the gravity always stay at right beneath the hips which makes body position in the chair always stable.



For those who have tonus and / or reflex movement in the chair, special attention by caregiver is needed in order to avoid falling down.



For those who have tendency of one-side-tilting, you may put Triangle Mat into the other side of the chair with careful attention.

## ■ LENGTH OF TIME -OF-USE



Even if it is so relaxing to use the Gravity Chair, using more than 30 minutes at the same position may cause a pain. We recommend you to change sitting position time to time and / or getting out of the chair within 30 minutes. For the people who have difficulty to change their body position by themselves, special attention by the caregiver is needed.



The Gravity Chair is mainly made of Urethane Foam and Polyurethane Leather Fabric. Ask at the nearest government office how to dispose the chair.



Use a firmly squeezed soft rag with mild soap to wipe off the covers to clean. DO NOT wash the cover with washing machine.

## How To Use Belt Holder Rubber Cord



Rubber cords on the edge of the arm rests are to hook the belt ends to make people easily sit into the chair without tangling the belts.

---

---

---

---



Pick up the seat belt and put the end of it through the rubber cord.

---

---

---

---



Make sure the seat belt is fully opened side ways so that the user can sit into the chair without being interrupted by the seat belt.

---

---

---

---