

Mannheimer

toilet and shower seat shell chair (MTS)

The Mannheim shower and toilet seat shell chair provides secure support and seating comfort. The spreading posture of the thighs, so typical when using the toilet, is encouraged. Slipping forward is prevented by a safety bar that can also be fitted with a hip retractor so that gentle pressure is applied to the lower abdomen - especially when the user leans forward. In addition, a backward tilt can be created with the help of the tilt adjustment. All this makes the seat pan chair a successful aid for toilet training.

The anatomical shape of the Ortho-Dynamic seat shell enables sitting without discomfort, even over longer periods of time. The comfort can be further increased if the seat is equipped with cushions. For patients, who require a seat according to the body for example, a milled seat made from closed-cell foam can be used. Four castors make it easy to push the seat shell chair over any toilet bowl. The equipment with a toilet bucket is also possible. The metal parts are made entirely of aluminum and stainless steel and the upholstery made of wet-strength foam. Thus chair can also be used as a shower chair.



Data	Size					
	2	3	4	5	6	7
Total weight	8 kg	9 kg	10 kg	12 kg	14 kg	16 kg
Total width	43	48	53	58	63	63
Total depth	65	65	65	70	70	70
Seat width rear	24	28	34	38	43	43
Seat width front	33	39	45	53	58	60
Width between arm rests	45	50	55	60	65	65
Seat height min.	48	48	48	48	48	48
Seat height max.	55	55	55	55	55	55
Seat depth	24	30	36	38	40	52
Back height	37	46	53	60	62	77
Toilet opening width (1)	12x17	14x20	16x23	18x26	20x29	20x29
Seat tilt forward (degrees)	0	0	0	0	0	0
Seat tilt back (degrees)	37°	37°	37°	37°	37°	37°

(1) Individual customizing possible



Size 4 and 6 compared

Basic configuration:

- Seat shell made of reinforced plastic
- Tilt adjustment (from size 3 with gas spring)
- Height adjustable
- Mobile on lockable castors, ø Größe 2 und 3: 75 mm, ø Größe 4-7: 100 mm
- Push handle
- headrest mount
- Maximal user weight: 120 kg

Accessories

- Headrest shell form
- Headrest with lateral support
- Headrest according to individual measures
- safety bar, padded
- hip retraction retainer
- hip belt
- Armrests, padded, height adjustable
- Abduction block, removable as splash guard
- Footrest, height adjustable
- Cutout reducer made of closed-cell foam
- Cushions for seat and back, cover in
- Toilet facility

Colours:

Seat shell RAL 9010 pure white (standard) other colours on request
underframe RAL 9006

The MTS is also available as a pure toilet attachment without a base. It can be equipped with almost any accessory.

Equipment options:



Headrest, push handle and hand release lever for tilting



Safety bar, seat and back cushions



Cutout reducer and back cushion with side support



Toilet facility



Foot release for tilting



Footrest height adjustable and foldable



MTS canted



Abduction block, removable as splash guard