

Multiroller



Buggy, Jogger, Offroader and much more

The Multiroller is designed for children with mild, moderate but also severe disabilities. These can be cerebral related disorders, muscle weaknesses or other handicaps.

Easy to lift

The basic set of the Multiroller with a sling seat weights less than 9 kg.

Advantage: It is often women who have to push the, lift it into the car or into the house. The low weight makes this work easier.

Easy to move

The Multiroller is designed as a high-tech product with a clever folding mechanism. The aim of the development was above all lightness, stability and German quality. The frame is particularly hard and rigid due to its special aluminium alloy. This means that the kinetic energy is converted into forward motion to a high and is not dissipated by twisting and rolling movements.

Advantage: Here too, it is mostly women who carry the larger and heavier child, even over uneven stretches. And then light rolling behaviour is a piece of additional Freedom aufgezehrt.

Universally applicable

The Multiroller can be used as a lightweight second supply, consisting only of the frame and a sling seat. With the adjustable backrest angle, it is a few degrees more sophisticated.

Numerous accessories for positioning and fixation of the child make it suitable more even for difficult clinical pictures. The special charm of the product is that it can be combined with other seating systems:



Buggy version with sling seat

- Individual seat shells
- Therapeutic seating systems from other manufacturers
- Many car seats of different brands and especially of RehaNorm (Carrot 3, Tony Plus).

Advantage: Of the numerous advantages only one example should be mentioned here: For a vacation trip, only a Multiroller and a Carrot 3 car seat are enough to transport the child.

Space-saving

As a buggy, jogger or off-roader, the Multiroller can be folded easily and requires little storage space, especially when the rear and front wheels are removed.

Good design

With this multitude of above-mentioned advantages, it is not a disadvantage if the product also looks good.

Conclusion

It is a buggy, it is a jogger, it is an offroader, it is a bicycle trailer and much more. It grows with your child and can carry up to 75 kg.

Colour

Silver metallic

Two frames - three front sets at the touch of a button



Offroader frame with four non-steering 16" wheels for forest and field and sand and beach



Buggy frame with two 16" rear wheels, two 10" swivel front wheels for use indoors or on flat ground.



Jogger frame with a non-steering 12" front wheel for sporty purposes such as skating or jogging. With a bicycle connection it becomes a bicycle trailer (see accessories)



Size 2 (as above but 5 cm higher)



Size 3 (as above but 5 cm wider)



Size 4 (as above but 5 cm higher and 5 cm wider)

Can be combined...



With jogger set and bicycle connection



Seat shell adaptor

...with many seats



With sling seat, footrest 1 and buggy set



With seat frame, footrest and offroad set



With seat shell, footrest and buggy set



With car seat Carrot 3, footrest and jogger set



With car seat Tony Plus, footrest 2 and buggy set



Also with corset



Seat frame...



...in many positions



And folded the entire Multiroller is very small

...with many accessories



Example: sun/rain canopy



Example: cargo bag



Example: head rest

Accessories

For solution 1: Sling Seat

- Sling seat made of spacer fabric blue or black
- Seat extension from 30 to 35 cm
- Back extension by 15 cm
- Headrest velcro
- Foot rest 1 fixed to the frame

For solution 2: Fixed seat

- Seat frame, with seat and back angle adjustable,
- seat depth 30 to 35 cm, cover: spacer fabric blue or black
- Back extension 15 cm
- Headrest velcro
- Footrest 2 fixed under the seat of the frame

For solution 3: Seat shell

- Plate for mounting seat shell adapters

For solution: Car seat

- Adapter für unterschiedliche
- Therapiesitze (auf Anfrage)

Für Lösung 5: Autositz

- Adaptor for Carrot car seat
- Adaptor for Tony Plus car seat
- Adaptors for other car seats (on request)

Other Accessories

- Additional buggy wheel set
- Additional jogger wheel set
- Additional offroad wheel set
- Bicycle connection (jogger wheel set necessary)



- Cushion for small children
- Trunk support with back plate
- Pommel on plate
- Seat wedge made to measure
- Harness
- 4- and 5-point-belt
- 4-point-shoulder-body-belt
- 2- and 4- point hip belt
- Groin strap
- Footrest straps
- Sun/rain canopy
- Leg/foot protection blanket
- Cargo bag under the buggy

Technical date	Narrow frame (dimensions in cm)		Wide frame (dimensions in cm)	
	Gr. 1 Drive wheels 16"	Gr. 2 Drive wheels 20"	Gr. 3 Drive wheels 16"	Gr. 4 Drive wheels 20"
Height total	102	107	102	107
Width total	68	68	73	68
Length total:				
- Buggy	116 bis 135	116 bis 135	116 bis 135	116 bis 135
- Jogger	136 bis 151	155 bis 170	136 bis 151	155 bis 170
- Offroad	128 bis 143	120 bis 135	128 bis 143	120 bis 135
Minimum turning circle:				
- Buggy	122	122	132	132
- Jogger	142	142	152	152
- Offroad	134	134	154	154
Height push handle	86 bis 116	86 bis 120	86 bis 116	86 bis 120
Footrest 1 W x L	38 x 23	38 x 23	43 x 26	43 x 26
Measures folded H x W x L	34 x 56 x 84	34 x 56 x 84	30 x 63 x 88	30 x 63 x 88
Maximal user weight	75 kg	75 kg	75 kg	75 kg
Weight with sling seat starting	9 kg	9 kg	10 kg	10 kg
Sling seat	58	58	58	58
- Height of back	33	33	33	33
- Depth of seat	30 – 38 (in front)	30 – 38 (in front)	35 – 43 (in front)	35 – 43 (in front)
- Width of seat	30 – 40 (at the top)	30 – 40 (at the top)	30 – (at the top)	30 – (at the top)
- Width of the back				
Integrated seat frame				
- Height of back	62	62		
- Depth of seat	30 - 35	30 - 35		
- Width of seat	38	38		
- Width of the back	38	38		
- Back adjustment in degrees	185° bis 225°	185° bis 225°		
- Seat adjustment in degrees	95° bis 135°	95° bis 135°		
- Lower leg length	20 – 30	20 – 30		
- Footrest 2 W x L	30 x 20	30 x 20		
- Footrest 2 adjustability in degrees	90° - 30°	90° - 30°		

- Persons with impaired vision can access this brochure on the website www.rehanorm.de and print it out in enlarged letters, if necessary, have it read aloud.
- Maintenance instructions, cleaning instructions and other technical specifications are contained in the user manual, which is available on the website www.rehanorm.de.
- When transporting in a vehicle, the child must not sit in the buggy but in a car safety seat. Also make sure that the Multiroller is safely stowed in the vehicle.